



### ScanPoint® Online

#### For BladderScan® Bladder Volume Instruments and the FloPoint® Elite Uroflow System

ScanPoint® Online offers expanded capabilities for BladderScan® BVI 6000 and BVI 9000 Series instruments and for the FloPoint® Elite Uroflow System.

Using ScanPoint® Online you can:

- View ultrasound measurements and images from patient exams
- Print exam results for patient records and for reimbursement (if applicable)
- Save and store images online and use the ScanPoint® automatic file backup service
- Ensure the optimum accuracy of your instrument by calibrating online, conveniently, in a matter of minutes and as often as you wish
- Enjoy software upgrades automatically



### ScanPoint® Local Client

#### For BladderScan® BVI 6100, BVI 6400 and BVI 9400 Bladder Volume Instruments and the FloPoint® Elite Uroflow System

If you are not ready or wired to take advantage of ScanPoint® Online bladder imaging, consider **ScanPoint® Local Client** (LC). With ScanPoint® LC:

- You can view, print and save images on your desktop.
- Calibration is available once per year via our loaner program.





### ScanPoint® Total Reliability<sup>SM</sup> Plan

The ScanPoint® Total Reliability Plan is included with your ScanPoint® subscription and provides :

- ScanPoint® online calibration for your BladderScan® BVI 6000 or BVI 9000 series instruments or annual factory calibration for your BVI 3000 instrument
- Technology and software assurance
- Instrument insurance
- Provision of a free loaner under our Next-Day Loaner<sup>SM</sup> program if your instrument needs to be evaluated
- Coverage of all shipping charges

While under continuous coverage from the ScanPoint® program, your instrument is warranted perpetually against outdated technology, which means we'll send you a new unit without charge when your instrument model is no longer manufactured and cannot be repaired. You'll also receive software updates for your BladderScan® instrument, via ScanPoint® Online or at no charge at an authorized service center.

Noninvasive, Accurate, Reliable, & Easy to Use

Benefit	ScanPoint® Online	ScanPoint® LC (Local Client)
		
<p>Choose the ScanPoint® service that best suits your practice</p>	<p>For the practice that:</p> <ul style="list-style-type: none"> <li>• Prefers the convenience and speed of Internet-based imaging services</li> <li>• Prefers the less hardware-intensive approach of online exam record storage</li> <li>• Prefers the convenience of calibration online and the immediacy of online technical support</li> <li>• Uses the BVM 6500 to measure ultrasound-estimated bladder weight (UEBW)</li> </ul>	<p>For the practice that:</p> <ul style="list-style-type: none"> <li>• Prefers access to exam records using their patient data management program, systems and reporting templates</li> <li>• May not be ready or wired to take advantage of ScanPoint® Online imaging services</li> <li>• Doesn't need the immediacy of online software updates or calibration. Can accommodate shipping the unit to Verathon® for calibration service</li> <li>• Doesn't require imaging services for the BVM 6500 instrument</li> </ul>
<p>Security and data backup</p>	<p>Automatically offers patient privacy protections (HIPAA compliant) Data backup is routine and automatically provided</p>	<p>Requires customer capability to ensure security and integrity of exam record storage Requires routine data backup by customer</p>
<p>View images</p>	<p>BVI 6100, BVI 6400, BVI 6500; BVI 9000 Series; FloPoint® Elite</p>	<p>BVI 6100, 6400; BVI 9400; FloPoint® Elite</p>
<p>Print images</p>	<p>BVI 6100, BVI 6400, BVI 6500; BVI 9000 Series; FloPoint® Elite</p>	<p>BVI 6100, 6400; BVI 9400; FloPoint® Elite</p>
<p>Save images</p>	<p>BVI 6100, BVI 6400, BVI 6500; BVI 9000 Series; FloPoint® Elite</p>	<p>BVI 6100, 6400; BVI 9400; FloPoint® Elite</p>
<p>Storage of data (exam records)</p>	<p>Via Internet Verathon® provides offsite storage at no additional charge</p>	<p>Via Local Client in local format Customer must make additional storage requirements</p>
<p>Calibration and technical support</p>	<p>Calibrate as often as desired via the Internet Verathon® can troubleshoot issues real time</p>	<p>Available once per year; a loaner is provided</p>
<p>Software</p>	<p>"Thin client" – No additional software required</p>	<p>Install CD provided</p>
<p>ScanPoint® upgrades</p>	<p>Easily accessible online</p>	<p>Available via mailed CD</p>
<p>Instrument upgrades</p>	<p>Device software updates available via ScanPoint® Online</p>	<p>Handled when the unit is sent to Verathon® for calibration</p>

US Patent No. 6569097, 6884217. INTL Patents Pending. BladderScan®, FloPoint®, ScanPoint®, Verathon® and Verathon Medical® are either registered trademarks or trademarks of Verathon Inc. in the USA and/or other countries. Copyright ©2008. All rights reserved.



**Corporate Headquarters**

2001 North Creek Parkway / Bothell, WA 98011-8218 USA  
425.867.1348 / 800.331.2313 / 425.883.2896 Fax

[www.verathon.com](http://www.verathon.com)