




Designed to Work
with GlideScope
Video Laryngoscopes

GlideRite Stylets—Enabling Quick Intubations



Designed to Work Together

GlideRite® Stylets are designed to complement angulated GlideScope® video laryngoscopes and enable quick placement of ET Tubes.

- Rigid design enhances Endotracheal Tube control for fast intubation
- Preformed to match the shape of GlideScope angulated blades
- Reusable and single-use options
- Options available to cover a broad range of patients and procedures
- Available in packs of 10

Single-use Stylets

Packaged for immediate use.



GlideRite Single-use Stylet – Small

Preformed, rigid biocompatible plastic

- Designed for use in ET Tubes I.D. 3.0 to 4.0
- Adjustable stopper



GlideRite Single-use Stylet – Medium

Preformed, rigid stainless steel

- Designed for use in ET Tubes I.D. 4.5 to 5.5
- Easy to use handle design for simple ET Tube securement and stylet removal



GlideRite Single-use Stylet – Large


Preformed, rigid stainless steel

- Designed for use in ET Tubes I.D. 6.0 and larger
- Easy to use handle design for simple ET Tube securement and stylet removal

Indicates single-use

Reusable Stylet

Compatible with high-level disinfectants and autoclave.



GlideRite Rigid Stylet

Preformed, rigid stainless steel

- Designed for use in ET Tubes I.D. 6.0 and larger
- Easy to use handle design for simple ET Tube securement and stylet removal

Schedule a demo with your Verathon representative or visit verathon.com/stylets

verathon

verathon.com

Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA

Tel: +1 800 331 2313
(USA & Canada only)

Tel: +1 425 867 1348
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.
Willem Fenengastraat 13
1096 BL Amsterdam
The Netherlands

Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited
Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia

Tel: 1800 613 603 (Within Australia)
Fax: 1800 657 970

Tel: +61 2 9431 2000 (International)
Fax: +61 2 9475 1201