




# Trusted Accuracy Meets Tested Durability



**BladderScan i10™**  
verathon



## BladderScan i10— The Next Generation Bladder Scanner

BladderScan i10™, powered by ImageSense® technology, delivers trusted accuracy for every user, accelerated workflow, tough and ready dependability, and certainty in your investment.

A partnership with Verathon® goes beyond the device itself. We deliver a comprehensive range of services tailored to your needs, giving you the power to maximize BladderScan® usage, minimize your cost of ownership, and improve patient care.




# Trusted Accuracy for Every User

Powered by ImageSense®, Verathon's groundbreaking artificial intelligence technology, BladderScan i10™ harnesses real-world data to consistently detect and define bladder location, size and shape across all pediatric and adult patient types and anatomies.

**± 7.5% on volumes greater than 100 mL<sup>1</sup>**

± 7.5 mL on volumes less than 100 mL



**CAUTI is largely preventable. Effective use of bladder scanner technology reduces unnecessary catheterization by nearly 50%<sup>2</sup>, lowering the risk of infection.**

**Successful CAUTI reduction initiatives depend on broad use and adoption by patient care staff of reliable and dependable bladder scanner technology.**

1. Accuracy specifications assume BladderScan i10 is being used in accordance with stated instructions, scanning a tissue-equivalent phantom.

2. Clinical Care Improvement Strategies, Preventing Catheter Associated Urinary Tract Infections, 2013 The Joint Commission.

# Accelerated Workflow

ImageSense®, combined with VMODE® probe technology, automatically measures the bladder volume within seconds.


- Large 10.1" (25.6 cm) optically bonded touchscreen reduces glare
- Interface uses intuitive graphics, automatic bladder location detection and BladderTraq™ contour tracking to guide the user through every step
- Simple charting tool easily transfers exam results to any electronic healthcare record system

Optimized for roaming through hectic hospital corridors, BladderScan i10 pairs with a quiet glide workstation.

- Single-handed probe docking keeps users focused on the patient – not on managing the equipment
- 360-degree handle enables easy positioning
- With both onboard charging and user-swappable battery, the system is always ready whenever and wherever you need it

## Simple EHR Charting

Scan the touchscreen directly, or scan the printout.



- Large display with intuitive graphics
- Tilt and height adjustable
- VESA quick mounting
- Drip resistant IPX2 sealing

- Rugged 3D probe
- Fully-sealed IPX4 protection
- Light and flexible cable
- Single-handed docking

360-degree handle for easy transport and positioning

- Aluminum column
- No-slip track mounting

Printer (optional)

Generous storage space for supplies

- AC adapter for onboard charging
- Cable management and protection

Low-profile base with quiet glide wheels

# Tough and Ready Dependability

BladderScan i10™ is designed and built for continuous operation in taxing hospital environments.


- Durable construction consisting of medical grade materials, metal gears and industrial-strength locking cables ensure reliable operation and high service availability
- No routine calibration is needed, even if you swap probes
- CaliScan™ probe diagnostics confirm probe function in the event of a drop
- BladderScan i10 passes a probe drop test that meets or exceeds a 4' (1.22m) drop to a steel surface using the MIL-STD 810H standard

BladderScan i10 is designed to meet the new realities of infection control.


- Custom-blended medical grade materials withstand fast-kill disinfectants
- Rugged, fully-sealed IPX4 probe and drip resistant IPX2 console ensure your precision electronics are protected

## Clean with Ease

Chemically resistant plastics used at each critical touchpoint.



# Product Highlights



**Accuracy for Every User**



**Accelerated Workflow**



**Tough and Ready Dependability**



**Flexible Power Management**



Click or scan to view 360°

# Certainty in Your Investment

You can keep out-of-pocket costs to a minimum over the life of your device, thanks to a comprehensive factory warranty and exceptional support.

Complimentary on-site training, software updates and fleet management services optimize usage. Technical support and customer service teams are on call to augment your on-site support.

Verathon® enhances the value of BladderScan i10™ with exceptional support, so you can have peace of mind knowing your investment is working every day to improve patient care.



## Comprehensive 5-Year Warranty

- 5-Year Factory Warranty covers all damage and failures associated with the probe and console
- Same-day loaner shipping when your order is processed by 2:00 PM PST
- “One call does it all”. If you need assistance, our experienced US-based customer service and full-line technical support teams are ready to help.
- Typical service turnaround within 5-days from receipt of instrument

## Optimized Uptime

Your customer experience is our top priority. Here are a few examples of how we’re making your life easier:

- No routine calibration needed
- Onboard diagnostics to support your preventative maintenance program
- Complimentary software upgrades allow you to update your instrument over time
- Quick transfer of configuration settings: Clone preferred device settings across the fleet
- Device configuration lockdown functionality ensures compliance with individual hospital regulations and best practices
- Service Partner Programs



**“Verathon’s on-site support has reinvigorated our bladder care program and helped us improve our CAUTI prevention practices.”**

Jonathan Kuehn, BSN, RN-BC  
Aspirus Wausau Hospital

# Training and Education

It takes more than a great device to improve patient care. You also need the know-how that enables you to fully integrate the device into your practice. Verathon® provides the education and training to enable clinical adoption of your BladderScan® purchase:

- On-site training
- Comprehensive training videos to educate on the product and bladder scanning skills
- In-service resources: posters, sign-in sheets, staff competency certification, and comprehensive overviews and quizzes

Let our team help you develop a plan that maximizes your BladderScan usage and streamline your operations so you can get the most out of your system.




# Accessories



Workstation  
0800-0631



Barcode Scanner  
0810-0291



Calibration Tank  
0620-1083



Tissue Phantom  
0620-0068



Printer  
0800-0640




Thermal Printer Paper (x 5)  
0800-0319



Ultrasound Gel  
0800-0005



Smart Battery Charger  
0400-0157



Rechargeable Battery  
0400-0164



Ultrasound Gel  
0800-0643  
20g single use packette,  
Box of 100

Learn more about the full BladderScan® suite of products,  
clinical resources and In-service instructional videos at [bladderscan.com](https://bladderscan.com)

---

**verathon**

[verathon.com](https://verathon.com)



Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA

Tel: +1 800 331 2313  
(USA & Canada only)

Tel: +1 425 867 1348  
Fax: +1 425 883 2896



Verathon Medical (Europe) B.V.  
Willem Fenengastraat 13  
1096 BL Amsterdam  
The Netherlands

Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia

Tel: 1800 613 603 (Within Australia)  
Fax: 1800 657 970

Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201